



MEDICAL

i3 trackApp

App for the digital filter change



INTENDED USE & FIELDS OF APPLICATION

The i3 trackApp is the ideal software solution for a digital, efficient, and, above all, effortless filter change. It enables to quickly plan, execute, and document the filter change in a structured and efficient way – providing a high degree of certainty and reliability. Data regarding the directions to the filters, the number of filters required, and the replacement dates are captured in full and are shown as a clear overview. The user-friendly handling by means of a smartphone or tablet enables the complete documentation of all water points.

FEATURES

Ensure a secure, simplified, and paperless documentation process

Digital planning on smartphone or tablet – even without an internet connection

Software compatibility with the operating systems Android and iOS

Simple and intuitive software application

Data is protected, secure, and stored digitally

Automatic directions and reporting

Safe personnel resources in the clinical environment





MEDICAL

TECHNICAL SPECIFICATIONS & REQUIREMENTS



APPLICATIONS

Filter change



PROVIDER

i3 Membrane GmbH



APP DATA SECURITY

Developed and hosted in Germany



COMMISSIONING

With the login data from i3 Membrane



REPORT

Transmission as a PDF file via email



APP FORMAT & COMPATIBILITY

Web app – compatible with all operating systems and devices (PC, smartphone, tablet), and features standard browsers (exception: Internet Explorer)



LANGUAGES

German
English



SCAN FUNCTION

Access to a functional camera is required

REGISTRATION & INSTALLATION



www.i3trackapp.de/en/registration

www.i3membrane.com/support

COMPATIBILITY

FILTER



i3 ONE



i3 TWO_{direct} / i3 TWO_{connect}



i3 FOUR



i3 FIVE

