



Catalog of Couplings

Your Companion
Life Science Separation






Table of contents

Couplings i3 ONE	pp. 4-11
Adaptors i3 ONE	pp. 12-13
Couplings i3 TWO	pp. 14-15
Couplings i3 FOUR	pp. 16-19
Couplings i3 FIVE	pp. 20-23
Overview of couplings	pp. 24-27

3

ONE

Couplings i3 ONE






Couplings i3 ONE

Water tap

Designation:	M22 F Connector ONE (AS)
Item number:	2060000
Description:	Quick-release coupling with M22x1 female thread AquaStopp function* included
Compatibility:	i3 ONE Point-of-Use water filter
Application:	Water tap with M22x1 male thread
Material:	Stainless steel (main body)
Scope of delivery:	1x M22 F Connector ONE (AS) 1x M22 flat gasket
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.




*The AquaStopp function prevents unintentional water leakage when a filter is not installed.



Couplings i3 ONE

Water tap

Designation:	M22 F TEMP Connector ONE (nAS)
Item number:	S206026
Description:	Quick-release coupling with M22x1 female thread for temporary use (max 150 days)
Compatibility:	i3 ONE Point-of-Use water filter
Application:	Water tap with M22x1 male thread
Material:	Polyethylene (PE) (main body)
Scope of delivery:	1x M24 M TEMP Connector ONE (nAS) 1x M22 F Adaptor 1x flat gasket
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.






Couplings i3 ONE

Water tap

Designation:	M24 M Connector ONE (AS)
Item number:	2060001
Description:	Quick-release coupling with M24x1 male thread AquaStopp function* included
Compatibility:	i3 ONE Point-of-Use water filter
Application:	Water tap with M24x1 female thread
Material:	Chrome plated brass (main body)
Scope of delivery:	1x M24 M Connector ONE (AS) 1x M24 flat gasket
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.




*The AquaStopp function prevents unintentional water leakage when a filter is not installed.



Couplings i3 ONE

Water tap

Designation:	M24 M TEMP Connector ONE (nAS)
Item number:	S206025
Description:	Quick-release coupling with M24x1 male thread for temporary use (max 150 days)
Compatibility:	i3 ONE Point-of-Use water filter
Application:	Water tap with M24x1 female thread female thread
Material:	Polyethylene (PE) (main body)
Scope of delivery:	1x M24 M TEMP Connector ONE (nAS) 1x M24 flat gasket
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.






Couplings i3 ONE

Water tap

Designation:	15/16"-27 M Connector ONE (nAS)
Item number:	2060009
Description:	Quick-release coupling with 15/16"-27 male thread AquaStopp function* not included
Compatibility:	i3 ONE Point-of-Use water filter
Application:	Water tap with 15/16"-27 female thread
Material:	Chrome plated brass (main body)
Scope of delivery:	1x 15/16"-27 M Connector ONE (nAS) 1x 15/16"-27 flat gasket
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.






Couplings i3 ONE

Water tap

Designation:	55/64"-27 F Connector ONE (nAS)
Item number:	2060019
Description:	Quick-release coupling with 55/64"-27 female thread AquaStopp function* not included
Compatibility:	i3 ONE Point-of-Use water filter
Application:	Water tap with 55/64"-27 male thread
Material:	Stainless steel (main body)
Scope of delivery:	1x 55/64"-27 F Connector ONE (nAS) 1x 55/64"-27 flat gasket
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.






Couplings i3 ONE

Water tap

Designation:	Sampling Adaptor
Item number:	S2060018
Description:	Single use adaptor for taking water samples from the quick-release coupling
Compatibility:	i3 ONE / i3 FOUR quick-release couplings
Application:	Water sampling
Material:	POM
Scope of delivery:	10x Sampling adaptor



Adaptors i3 ONE






Adaptors i3 ONE


Water tap

Item number	Specifications*
GX332310001	10PK-9/16"-24 MALE X 55/64"-27 MALE
GX333310001	10PK-11/16"-24 MALE X 55/64"-27 MALE
GX336310001	10PK-15/16"-27 MALE X 55/64"-27 MALE
GX337310001	10PK-3/8" MALE X 55/64"-27 MALE
GX339310001	10PK-13/16"-27 MALE X 55/64"-27 MALE
GX341310001	10PK-3/4"-27 FEMALE X 55/64"-27 MALE
GX344310001	10PK-1/16"-27 MALE X 55/64"-27 MALE
GX349310001	10PK-13/16"-24 FEMALE X 55/64"-27 MALE
GX362310001	10PK-3/8"-18 NPT FEMALE X 55/64"-27 MALE
GX365310001	10PK-1/4"-18 NPS MALE X 55/64"-27 MALE
GX367310001	10PK-5/8"-24 MALE X 55/64"-27 MALE
GX377310001	10PK-5/8"-27 MALE X 55/64"-27 MALE

*Available only in packs of ten.



Couplings i3 TWO






Couplings i3 TWO


Shower

Designation:	G1/2" M Connector TWO (AS)
Item number:	2060002
Description:	Quick-release coupling with G1/2" male thread AquaStopp function* included
Compatibility:	i3 TWOconnect POU water filter
Application:	Shower hose with G1/2" female thread
Material:	PSU (main body)
Scope of delivery:	1x G1/2" M Connector TWO (AS) 1x 1/2" flat gasket
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.



*The AquaStopp function prevents unintentional water leakage when a filter is not installed.

Couplings i3 FOUR






Couplings i3 FOUR

Inline

Designation:	G1/2" M Connector FOUR (nAS)
Part number:	2060060
Description:	Quick-release coupling with G1/2" female thread AquaStopp function* not included
Compatibility:	i3 FOUR Inline filter
Application:	Cleaning guns with G1/2" female thread
Material:	Stainless steel
Scope of delivery:	1x G1/2" M Connector FOUR (nAS) 1x G1/2" flat gasket
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.




*The AquaStopp function prevents unintentional water leakage when a filter is not installed.



Couplings i3 FOUR

Inline

Designation:	G1/2" F Connector FOUR (AS)
Item number:	2060061
Description:	Quick-release coupling with G1/2" female thread AquaStopp function* included
Compatibility:	i3 FOUR Inline filter
Application:	Cleaning guns with G1/2" male thread
Material:	Stainless steel
Scope of delivery:	1x G1/2" F Connector FOUR (AS) 1x G1/2" flat gasket
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.



*The AquaStopp function prevents unintentional water leakage when a filter is not installed.



Couplings i3 FOUR Inline

Designation:	Adaptor Lock
Item number:	2060071-10
Description:	Safety device against unintentional activation of quick-release coupling
Compatibility:	i3 FOUR Inline filter
Application:	Securing inline quick-release couplings
Material:	ABS
Scope of delivery:	10x adaptor lock




open



closed

Couplings i3 FIVE







Couplings i3 FIVE

Inline

Designation:	G1/2" M Connector FIVE (AS)
Item number:	2060050
Description:	Quick-release coupling with G1/2" male thread AquaStopp function* included
Compatibility:	i3 FIVE Inline filter
Application:	Endoscopy AER; birthing tubs; flex hoses with G1/2" female thread
Material:	POM (main body)
Scope of delivery:	1x G1/2" M Connector FIVE (AS)
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.




*The AquaStopp function prevents unintentional water leakage when a filter is not installed.




Couplings i3 FIVE

Inline

Designation:	G1/2" F Connector FIVE (AS)
Item number:	2060054
Description:	Quick-release coupling with G1/2" female thread AquaStopp function* included
Compatibility:	i3 FIVE Inline filter
Application:	Endoscopy AER; birthing tubs; flex hoses with G1/2" male thread
Material:	POM (main body)
Scope of delivery:	1x G1/2" F Connector FIVE (AS)
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.




*The AquaStopp function prevents unintentional water leakage when a filter is not installed.



Couplings i3 FIVE

Inline

Designation:	Wall Mount FIVE
Item number:	2060052
Description:	Wall mount for the i3 FIVE
Compatibility:	i3 FIVE Inline filter
Application:	Bracket for mounting i3 FIVE inline filter
Material:	Santopren (main body)
Scope of delivery:	1x Wall Mount FIVE 1x mounting kit
Installation:	Assembly and installation of the components must be carried out by qualified technical personnel.



Overview of Couplings i3 ONE




i3 ONE Water tap		
Image	Designation	Item number
	M22 F Connector ONE (AS)	2060000
	M22 F TEMP Connector ONE (nAS)	S206026
	M24 M Connector ONE (AS)	2060001
	M24 M TEMP Connector ONE (nAS)	S206025
	15/16"-27 M Connector ONE (nAS)	2060009
	55/64"-27 F Connector ONE (nAS)	2060019
	Sampling Adaptor	S2060018

Overview of Couplings

i3 TWO



i3 TWO Shower		
Image	Designation	Item number
	G1/2" M Connector TWO (AS)	2060002



Overview of Couplings

i3 FOUR



i3 FOUR Inline		
Image	Designation	Item number
	G1/2" M Connector FOUR (nAS)	2060060
	G1/2" F Connector FOUR (AS)	2060061
	Adaptor Lock	2060071-10

Overview of Couplings

i3 FIVE



i3 FIVE Inline		
Image	Designation	Item number
	G1/2" M Connector FIVE (AS)	2060050
	G1/2" F Connector FIVE (AS)	2060054
	Wall Mount FIVE	2060052



i3 Membrane Corporation
13801 Reese Blvd West, Ste 310
Huntersville NC, 28078
USA

i3 Membrane GmbH
Christoph-Seydel-Str. 1
DE-01454 Radeberg
Germany

i3 Membrane GmbH | Global
medical@i3membrane.de

i3 Membrane Corporate | U.S.
support@i3membrane.com

Visit our website at
www.i3membrane.com



Your Companion
Life Science Separation

